

## Composition

### P07-03110 ITS solution I (100x)

	Components	mg/L
Inorganic Salts	Calcium chloride anhydrous	200.00
	Magnesium sulfate x 7H <sub>2</sub> O	200.00
	Potassium chloride	400.00
	Sodium chloride	6,800.00
	Sodium dihydrogen phosphate x H <sub>2</sub> O	140.00
	Sodium selenite	0.500
Other Components	D(+)-Glucose anhydrous	1,000.00
	Insulin human rec.	500.00
	apo-Transferrin human	500.00